



Product Catalog



SIRIUS +

Topographer & Tomographer

Page: 3-14



MS-39

Anterior Segment OCT

Page: 15-31



ANTARES

Topographer

Page: 32-40



MODI-02

Topographer

Page: 41-44



COBRA^{HD}

Retinal Camera

Page: 45-50



COBRA^{PLUS}

Retinal Camera

Page: 51-54



OSIRIS-T

Aberrometer & Topographer

Page: 55-58



OSIRIS

Aberrometer

Page: 59-62



SL-9800

Slit Lamp

Page: 63-73



SL-9900

Slit Lamp

Page: 63-73



SL-9800 ZOOM

Slit Lamp

Page: 63-73



SL-9900 ZOOM

Slit Lamp

Page: 63-73



VL-Z8

Digital Camera Adaptor

Page: 70



Z-900

Tonometer

Page: 72



Z-800

Tonometer

Page: 73



SL-9900 ELITE

Slit Lamp

Page: 76-85



HR-CAMERA

Digital Video Camera

Page: 82-83



LIFT -01&02

Electrical Tables

Page: 86-91



POLARIS

Tear Film Analysis

Page: 92-94



PERSEUS

Endothelial Microscope

Page: 95-100



SIRIUS⁺ SCHEIMPFLUG

Tomographer & Topographers



MADE IN ITALY



OMOGRAPHER




TOPOGRAPHER

ABOUT SIRIUS+



SIRIUS examinations provides an accurate measurement of pupil diameter in scotopic, mesopic and photopic conditions. When combined with the corneal map they can be used for refractive surgery planning and follow up. All biometric measurements of the anterior chamber are calculated using 25 sections from the cornea.






PRODUCT OVERVIEW

INCLUDED:

1. Chinrest Module
2. Head Rest
3. Capturing Channel
4. Resolution (1280x960)
5. Chinrest Cup
6. Knob Adjuster
7. Patient's Handle
8. Device Blocking Knob
9. Capturing Trigger Joystick
10. Slide Guide Guards

LIFT

**Electrical Table
02**



OPTIONAL:

11. Tools Storage
12. HP Monitor
13. Adjustable Stand
14. Custom Colors Available
15. Single & Double Work Stations
16. Table Height Adjuster
17. Elevation Column
18. Cogged Wheels



LIFT

**Electrical Table
01**




SCHIEMPFLUG TOMOGRAPHY

100 HR
Corneal Sections





Corneal Thickness


SIRIUS is a fully featured multi-functional placido disk topographer and Scheimpflug's tomographer with a dedicated software designed to help in the detection and analysis of Dry Eyes.



Anterior, posterior and
total refractive power
map reports



Anterior chamber depth
and posterior corneal
surface elevations




Tangential and axial curvature
of anterior and posterior
corneal surface

PRECISE SURFACES

SIRIUS is able to obtain the accurate measurement of elevations, curvature, power and thickness for the whole cornea surfaces over a 12mm in diameter.



The most common uses are for refractive and cataract surgery, and IOL calculation module is available. Photos and videos with blue & white diffuse light.



ILLUMINATION

Placido Disk	LED @400-700nm
Scheimpflug	LED @475nm UV-free
Pupillography	LED @940nm
Auxiliary Lighting	LED 400-700nm
Fluorescence Lighting	LED @470nm



PHOENIX

SOFTWARE

CONNECT YOUR TWO DEVICES

SIRIUS+ uses the Phoenix Software platform allowing patient data to be saved for future review and analysis, shared by all CSO devices. Enables both comfort and flexibility with CSO's optional single and double electrical tables designed for all ophthalmic instruments.



Windows




Simk						n 1.3375	
Simk	53.72 D			R	6.28 mm		
Flat Simk	51.82 D	155°		R1	6.51 mm		
Steep Simk	55.62 D	65°		R2	6.07 mm		
Astig	3.81 D	65°		e ² (-Q)	1.86		
Posterior Axial Curvature							
Mean K	-8.36 D			R	4.78 mm		
Flat K	-7.97 D	155°		R1	5.02 mm		
Steep K	-8.75 D	65°		R2	4.57 mm		
Astig	-0.78 D	65°		e ² (-Q)	2.08		
Total Corneal Power IOL (Ray Traced)							
Mean TCPIOL	51.88 D			Central	52.45 D		
Flat TCPIOL	50.13 D	155°		Mid	48.36 D		
Steep TCPIOL	53.64 D	65°		Periph	45.64 D		
Astig	3.51 D	65°					
Pachymetry							
o Thinnest	460 μm			location x,y	0.64 mm	-0.48 mm	
Central	498 μm			CCT	477 μm		
Mid	581 μm						
Periph	671 μm						
Anterior Axial Curvature Zones							
Central	54.46 D			6.20 mm			
Mid	48.19 D			7.00 mm			
Periph	44.09 D			7.66 mm			
Kmax	60.54 D			5.57 mm	location x,y	-0.01 mm	-0.80 mm
Anterior Chamber and Biometry							
WTW, N-T	11.61 mm			Mean Angle	36.6 °		
ACV	119 mm ³			Kappa Dist	0.50 mm		
AQD	2.95 mm			ASL endo	n/a		
+ Pupil Diam	2.46 mm			location x,y	0.13 mm	0.48 mm	
Keratocornus Probability							
KPI	100.0 %			Kprob	100.0 %		
CLMIaa	9.31 D			PPK	100.0 %		




ANTERIOR SEGMENT DIAGNOSTICS

Anterior and posterior corneal topography information are available for diagnosis, for refractive/ cataract pre-operative planning or for follow up purposes.

Pupillography



Dry Eye Report




Tear meniscus & Ocular Redness



RGP lens fitting



Tear Film Analysis



Meibomian Glands Analysis



Densitometry



Topographic Maps & Summary



What does it do?




- Map of Anterior Chamber Depth & Analysis
- Refractive & Cataract Surgery Planning Tools
- IOL calculation with Ray Tracing Techniques
- Automatic Calculation of Iridocorneal Surface and works with already treated eyes (e.g. LASIK)
- Advanced Topography Ring Editing System
- Corneal Pachymetry (12 mm diameter)
- Contact Lenses Application Module
- OPD Analysis and Visus Simulation
- Scheimpflug's Images Comparison
- Summary of Acquisition Reliability
- Intra Stromal Rings Summary
- Phoenix Software Platform
- IOP Correction Formulas
- 4 maps + image summary
- 24 Rings Placido's Disk
- Contact Lenses Autofit
- Keratoconus Summary
- Corneal Aberrometry
- Glaucoma Summary
- Videokeratoscope
- Cataract Summary




DRY EYE REPORT OD

Thanks to the new color camera, allows the measure of tear film break-up time, meibomian glands analysis, conjunctival, limbar hyperemia and tear meniscus height. In addition, all functionalities merge together for all partial score to provide a complete Dry Eye report for comprehensive assessment of the patients corneal condition and helping the diagnosis of the Dry Eye Disease DED.

Reset  Compile 	
Questionnaire	OSDI Mild
Reset  Grade  Open  Acquire 	
Ocular redness	Conjunctival Grade 2
Reset  Open  Acquire 	
Meibomian glands	Glands loss 7,1%
Reset  Open  Acquire 	
Tear meniscus	TM height 0,35mm
Reset  Open  Acquire 	
NiBuT	NiBuT 13,0sec
Reset  Import 	
Osmolarity	275 Normal







ADVANCED ANALYSIS OF THE TEAR FILM

Placido disk technology allows for the advanced analysis of the tear film, such as NI-BUT (Non Invasive Break-up Time).


CORNEAL ABERROMETRY

Offers a complete overview of the corneal contribution to the vision. Anterior, posterior or total corneal aberration are selectable for several pupil diameters. The OPD/ WFE map and the simulated vision functions (Spot Diagram, PSF, MTF, Image convolution) help the clinician understanding and explaining the visual discomfort to the patient.



CONTACT LENSES APPLICATION MODULE

A contact lens fitting module is available which simulates the fit of rigid lenses based on an internal database of many lens manufacturers.





PUPILLOGRAPHY

Sirius has built-in pupillography with pupil measurement of scotopic (0.04 lux), mesopic (4 lux), photopic (50 lux) conditions and in dynamic mode. Knowing the center and the diameter of the pupil, is essential for many clinical procedures which seek to optimize vision quality.

MEIBOGRAPHY


Meibomian glands can be viewed under infrared light once the image is captured, you can use the software to aid in the analysis of the condition of the glands



IOL CALCULATION MODULE

(OPTIONAL)

This module is based on Ray-Tracing techniques, regardless of the state of the cornea (untreated or previously treated for refractive purposes), provides the calculation of the spherical and toric power of the intraocular lens.





INTRASTROMAL RINGS

On the basis of the pachymetry map and corneal alti - metric data, SIRIUS allows for intrastromal rings system planning, with variable options for the correction of refractive defects and some forms of keratoconus.

GLAUCOMA SUMMARY

Aberrometric analysis offers a complete overview of the corneal aberrations. It is possible to select the contribution cornea of the anterior, posterior or total for different pupil diameters. The OPD/WFE maps and the visual simulations (PSF & MTF) can help the clinician in understanding or explaining the patient's visual problems.



KERATOCONUS SCREENING

Keratoconus screening provides the clinician with important information about the patients cornea. This will help prevent complications associated with ectasia before corneal surgery is undertaken.



PRODUCT SPECIFICATIONS

MADE IN ITALY

Technical Data

Data Transfer	USB 3.0
Power Supply	External power source 24 VCC
Power net cable	IEC C14 plug
Dimensions (HxWxD)	515 x 315 x 255mm
Weight	7 Kg
Chin rest movement	70mm \pm 1mm
Base movement (xyz)	105 x 110 x 30mm
Working distance	74mm

Light Source

Placido Disk	LED @400-700nm
Scheimpflug	LED @475nm UV-free
Pupilligraphy	LED @940nm
Fluoresceine Lighting	LED @470nm
Auxiliary Lighting	LED 400-700nm

Topography

Placido rings	22
Topographic Coverage	12mm
Measured Points	Class A according to UNI EN ISO 19980-2012
Dioptric Measurement range	1D to 100DT
Measurement accuracy	Class "A" as per "ISO19980:2005 (E)
Power frequency	(50/60Hz) magnetic field IEC 61000-4-8
Power cable	Four-core cable conductors



Accessories

Light Diffuser Filter	Auxiliary illumination & magnetic
Yellow Barrier Filter	Magnetic 530 nm filter
Additional Lens	Magnetic -6D lens
Calibration Tool	R8 mm calibration tool



MS-39

Anterior Segment OCT



MADE IN ITALY



VL *Anterior Segment OCT* MS-39




The MS-39 is the most advanced corneal map device using a SD-OCT and Placido disk corneal topography to obtain measurements of the anterior segment of the eye.

- MS-39 provides topographic maps:
 - Tangential curvature (anterior and posterior)
 - Sagittal curvature (anterior and posterior)
 - Elevation (anterior and posterior)
 - Refractive power (equivalent, anterior and posterior)
 - Corneal thickness
 - Epithelial thickness
 - Anterior chamber depth

- In addition to anterior segment clinical diagnostics, MS-39 can be used in corneal surgery for refractive surgery planning.

- MS-39 is known to produce the highest resolution and clarity of the cross-sectional images, with a 16 mm diameter, along with the many details of the cornea structure and layers.




OCT - Behind the Eyes



MS-39 can measure and detect on the Placido image so that height, slope, and curvature data can be calculated using the arc-step method with conic curves. Profiles of the anterior cornea, posterior cornea, anterior lens, and iris are derived from the SD-OCT scans.

Data for the anterior surface from the Placido image and SD-OCT scans are merged using a proprietary method. All other measurements for internal structures (posterior cornea, anterior lens, and iris) are derived solely from SD-OCT data.





USB CONNECTION

The USB connection between the device and the PC enables a fast and easy transfer of the images.



AUTOMATIC MEASUREMENTS

MS-39 combines a Placido disk corneal topography, with high resolution OCT based anterior segment tomography.



HIGH PRECISION JOYSTICK

MS-39 captures footage with the advanced manual acquisition and electronically guided control joystick, which guarantees high precision and repeatability of all measurements.







Phoenix Software

MADE IN ITALY




The MS-39 works perfectly with the advanced PHOENIX software. This program enables comfortable working, by connecting all of your diagnostic instruments with a powerful patient database, giving you an extraordinarily effective work station.




EPITHELIAL MAP

Epithelial thickness mapping can help us identify early keratoconus while screening patients for corneal refractive surgery. Epithelial thickness mapping provides additional supporting evidence in the diagnosis, which allows us to offer alternative treatment options.






Thickness Map



Placido Imaging

The MS-39 Anterior Segment OCT includes the innovative and advanced measurement of the epithelial layer. The epithelial masking effect is known for its morphology which is very useful assess abnormalities of the corneal surface.





CRYSTALLINE BIOMETRY

Lenticular crystalline biometry on non cataractous lenses has been studied by means of Scheimpflug photography and digital image analysis.




Anterior Lens Surface



Crystalline Lens

In order to more accurately determine and refine the intra-ocular lens calculation, the device provides a capture mode to measure the crystalline lens thickness, its distance from the cornea and its equator.

The crystalline lens of the eye is the only human organ which undergoes a steady increase in size due to growth throughout the life of the individual.



ADVANCED ANALYSIS OF THE TEAR FILM

Placido disk technology allows for the advanced analysis of the tear film, such as NI-BUT (Non Invasive Break-up Time).





Placido Disk



Meniscus Evaluation


VIDEOKERATOSCOPY MODULE

- Tear film examination with white light
- Tear film examination with fluorescein
- Break-up time measurement.
- Examination of tear layers.
- Examination of rigid LAC clearance with fluorescein




PUPILLOGRAPHY

Sirius has built-in pupillography measurement software. Knowing the center and the diameter of the pupil, is essential for many clinical procedures which seek to optimize vision quality.




3D Scheimpflug




Contour Detection

PUPILLOMETRY MODULE


- Pupillometry with scotopic light to determine pupil maximum extension and optic zone diameters for treatment settings.
- Pupillometry with mesopic light (4 lux)
- Pupillometry with photopic light (50 lux)
- Dynamic pupillometry, starting with over 400 lux and switching off the light source so that the pupil can dilate to its maximum extension
- Evaluation of pupil decentralization from the corneal vertex and calculation of the pupil centre during dilation.




KERATOCONUS SCREENING



Keratoconous screening provides the clinician with important information about the patients cornea. This can help prevent complications associated with ectasia before corneal surgery is undertaken.



Klyce / Maeda Indices



Rabinowitz Values

Keratoconus is a non-inflammatory, progressive thinning process of the cornea. It is a relatively common disorder of unknown etiology that can involve each layer of the cornea and often leads to high myopia and astigmatism.

IOL CALCULATION MODULE


This module is based on Ray-Tracing techniques, regardless of the state of the cornea (untreated or previously treated for refractive purposes), provides the calculation of the spherical and toric power of the intraocular lens.



Post Refractive




Toric Results



MEASURED POINTS

31,232
(anterior surface)

25,600
(posterior surface)



INTRASTROMAL RINGS

On the basis of the pachymetry map and corneal altimetric data, SIRIUS allows for intrastromal rings system planning, with variable options for the correction of refractive defects and some forms of keratoconus.



Intracorneal Ring




Surgical Incision

PROS:

- For patients with progressive keratoconus,
- Reversible
- Patients may achieve a better fit and hence a more comfortable fit with their contact lenses


CONS:

- Postoperative discomfort/pain for 1-4 days
- Vision may not improve or change
- Fluctuating vision for weeks to months after the procedure



RESHAPED CORNEA





GLAUCOMA SCREENING

The MS-39 device allows the glaucoma screening and gives the measurement of irido-corneal angles AOD, TISA and corneal pachymetry.




Healthy



Glaucoma


These values are the most common IOP correction formulas, useful to diagnose the eyes optic nerve and intraocular pressures to detect any diseases which can be due to the conformation of the anterior chamber.






CORNEAL ABERROMETRY

Aberrometric analysis offers a complete overview of the corneal aberrations. It is possible to select the contribution of the anterior, posterior or total cornea for different pupil diameters. The OPD/WFE maps and the visual simulations (PSF & MTF) can help the clinician in understanding or explaining the patient's visual problems.



Wavefront Imaging



OPD Maps

Assisted manual acquisition advanced ring editing system available maps:

1. Sagittal curvature map
2. Angular curvature map
3. Altimetry
4. Refractive power
5. Gaussian curvature map

MS-39

PRODUCT OVERVIEW





- | | | |
|---------------------------------|-------------------------------|-----------------------------------|
| 1. Chin Rest Module | 7. Table Drawer | 13. Chin Rest Knob Adjuster |
| 2. Head Rest | 8. Chin Rest Support | 14. Joystick (Capturing Trigger) |
| 3. Shooting Channel | 9. Electrical Table Adjuster | 15. Power Supply Connector |
| 4. Chin Rest Cup | 10. Device Locking Knob | 16. Device Supply Cable |
| 5. Personal Computer (Optional) | 11. Instrument / Placido Disk | 17. Slide Guard Guards |
| 6. Patient's Handle | 12. Paper for Chin Rest | 18. LIFT 02 - Tabletop (Optional) |





- | | | |
|--|----------------------------------|---------------------------------|
| 1. Monitor (Optional) | 6. Data Plate | 14. One Column Table (Optional) |
| 2. USB Connection (Device + PC) | 7. Power Supplier | 15. Power Supply Connector |
| 3. Power Cable Connection
(Power Supplier + Device) | 8. Supply Control Light | 16. Power Sockets |
| 4. Power Cable Connection
(Electrical Table + Power Supplier) | 9. ON / OFF Switch | 17. Two Column Table (Optional) |
| 5. Power Cable Connection
(Electrical Table + Power Supply) | 10. Power Supply Out Connector | 18. Table Placement |
| | 11. Power Supply Mains Connector | |
| | 12. Power Supply Cable | |
| | 13. Power Supply Cable (Out) | |

Not Included (Optional)



Not Included (Optional)

Not Included (Optional)





Technical Data

Data Transfer	USB 3.0
Power Supply	External power source 24 VCC In: 100-240Vac 50/60Hz - 2A - Out: 24Vdc - 100W
Power net cable	IEC C14 plug
Dimensions (HxWxD)	505 x 315 x 251mm
Weight	10.4Kg
Chin rest movement	70mm \pm 1mm
Minimum height of the chin cup from the table	23cm
Base Movement (xyz)	105 x 110 x 30mm
Working distance	74mm

Light Sources

Placido disk illumination	Led @635nm
OCT source	Led @845nm
Pupillographic illumination	Led @635nm

Topography

Placido disk rings	22
Measured points	31232 (anterior surface) 25600 (posterior surface)
Topographic covering	10mm
Dioptric measurement range	1D to 100D
Measurement accuracy	Class A according to UNI EN ISO 19980-2012

Section

Image field	16mm x 8mm
Axial resolution	3.6 μ m (in tissue)
Transversal resolution	35 μ m (in air)
Image(s) resolution	Keratotomy (640x480) + 25 radial scans on a 16mm transversal field (1024 A-scan) - Section: on 16mm (1600 A-scan) on 8mm (800 A-scan)
Operating system	Windows 10 (64 bit)




ANTARES


Topographer




HIGH PRECISION ACQUISITION



DATA AND ANALYSIS



- Pupillography
- Fluorescein Imaging
- Lipid layer pattern imaging
- Advanced Analysis of the Tear Film
- Contact Lenses Application Module
- Tear meniscus height imaging and measurement
- Meibomian glands analysis, imaging & screening
- Keratoconus screening and detection
- High-resolution color video camera
- Advanced ring editing system
- 24 rings Placido's disk
- Videokeratoscope



CORNEAL MEASUREMENTS

ANTARES topography function provides information about the curvature, elevation and refractive power of the cornea. Also provides many parameters to aid in the diagnosis and monitoring the corneal surface.



PUPILLOMETRY MODULE

- Pupillometry with scotopic light to determine pupil maximum extension and optic zone diameters for treatment settings.
- Pupillometry with mesopic light (4 lux)
- Pupillometry with photopic light (50 lux)
- Dynamic pupillometry, starting with over 400 lux and switching off the light source so that the pupil can dilate to its maximum extension
- Evaluation of pupil decentralization from the corneal vertex and calculation of the pupil centre during dilation.

CORNEAL TOPOGRAPHY MODULE

- Assisted manual acquisition
- Advanced ring editing system
- Available maps:
 1. Sagittal curvature map
 2. Angular curvature map
 3. Altimetry
 4. Refractive power
 5. Gaussian curvature map



VIDEOKERATOSCOPY MODULE

- Tear film examination with white light
- Tear film examination with fluorescein
- Break-up time measurement.
- Examination of tear layers.
- Examination of rigid LAC clearance with fluorescein

LIFT 01

Electrical Table - Not Included



Optional Tabletops

Available Colors

PRODUCT DESCRIPTIONS

ANTARES provides a multi topography report from the images captured from the built-in digital video camera. This device includes an editing software which allows you to edit edge position for proper reconstruction on all distorted surfaces.



ANTARES

Phoenix Software

MADE IN ITALY





ANTARES works perfectly with the advanced **PHOENIX** software. This program enables comfortable working, by connecting all of your diagnostic instruments with a powerful patient database, giving you an extraordinarily effective work station.



KERATOCONUS SCREENING


Keratoconous screening provides the clinician with important information about the patients cornea. This can help prevent complications associated with ectasia before corneal surgery is undertaken.



CONTACT LENSES APPLICATION MODULE

A contact lens fitting module is available which simulates the fit of rigid lenses based on an internal database of many lens manufacturers.

Lens Module




Anterior Segment




PUPILLOGRAPHY

ANTARES has built-in pupillography measurement software. The measurement of the pupil in scotopic (0.04 lux), mesopic (4 lux), photopic (50 lux) conditions and in dynamic mode. Knowing the center and the diameter of the pupil, is essential for many clinical procedures which seek to optimize vision quality.



Infrared Pupillograph

Pupil Detection




MEIBOGRAPHY

Meibomian glands can be viewed under infrared light once the image is captured, you can use the software to aid in the analysis of the condition of the glands

Infrared Imaging

Transillumination





ADVANCED ANALYSIS OF THE TEAR FILM

Placido disk technology allows for the advanced analysis of the tear film, such as NI-BUT (Non Invasive Break-up Time).

Placido Disk

Meniscus Evaluation






DRY EYE REPORT

Based on the Ocular Surface Disease Index questionnaire (OSDI), limbal and conjunctival hyperaemia, Meibomian glands analysis, tear meniscus analysis, NIBUT, and tear osmolarity, calculated merging together all partial scores, provides an overall evaluation of the clinical condition of the patient for a comprehensive diagnosis of the dry eye disease.

ANTARES

DEVICE DIAGRAM



1. Capturing Channel
2. Chin Rest Module
3. Head Rest
4. Calibration Tool
5. Chin Rest Cup
6. Knob Adjuster
7. Patient's Handle
8. Device Blocking Knob
9. Instrument with Placido's Disk
10. Joystick with Capturing Trigger
11. Slide Guide Guards
12. Power Supply Cable
13. Power Supplier

Measurements

Operation Distance	74 mm from corneal vertex
Number of Rings	24
Number of Measuring Points	6144 (24x256)
Number of Points Analysed	Over 100000
Diameter of the corneal area covered	0.4 to over 9.6 mm of Diameter
Dioptries Measuring Arc	1 to 100 D
Size (HxWxD) mm	505x315x251 mm
Weight	6.2 kg

Applicable Lighting

Placido's LED lighting	White LED
Fluorescein LED lighting	Blue LED 460 nm
Pupillometry	LED lighting IR LED 875 nm

Notes

Accuracy and repeatability error	Class "A" as per "ISO19980:2005 (E)
Power supply	24V DC external power supply unit
Input power supply unit	90-264 V AC: - 47/63 Hz Max 0.9 A OUTPUT: 24 V DC - 2 A
Power frequency	(50/60Hz) magnetic field IEC 61000-4-8
Power cable	Four-core cable conductors
Computer connection	USB3 Type A cable



MODI-02

Topographer


LIFT 02

Electrical Table - Not Included



Optional Tabletops


Available Colors



MODI-02


Topographer

1. Instrument with Placido's Disk
2. Chinrest Module
3. Headrest
4. Capturing Channel
5. Chinrest
6. Joystick with Capturing Trigger Button
7. Base Slides Locking Knob
8. Firewire Port
9. Patient's Handle
10. Glide Slides Guards



Features:

- Ergonomic design, high quality optics and precision mechanical parts
- MODI is an electro medical system for the detection, capturing and digital processing of an image of the cornea.
- MODI is a device that allows "live" shooting on the computer monitor.
- Management and control software including cornea measurement.
- High-resolution monochromatic video camera
- Multiple maps comparison in a single display window
- Placido's disk with 24 rings



Product Introduction

MODI-02 Topographer




Corneal Topographer

MODI technology is featured through multi-functional corneal topographer. It has a dedicated software designed to help in the analysis of the Eye.




Software Phoenix

Phoenix Software is designed for taking fine and detailed digital images, and sorting files. Software run in all Windows XP, 7, 8, and 10. (Included)




Pupillography

Antares measures the pupil in scotopic (0.04 lux), mesopic (4lux), photopic (50 lux) conditions and in dynamic modality is fast and simple.




Videokeratoscope

Antares has a white light to capture color images and videos as well as cobalt light blue light for analysis of contact lenses fitting with fluorescein.



Keratoconous Screening

The LED technology has an estimated life of 5000 of continuous use. Image colour temperature is maintained at all illumination levels.



Contact Lenses Module

The MODI application contact fitting module stimulates the fit of rigid contact lenses based on an internal database lens manufacturers.

MODI-02

Topographer



Product Specifications

Measurements

Operation Distance	74 mm from corneal vertex
Number of Rings	24
Number of Measuring Points	6144 (24x256)
Number of Points Analysed	Over 100000
Diameter of the corneal area covered	0.4 to over 9.6 mm of Diameter
Dioptres Measuring Arc	1 to 100 D
Size (HxWxD) mm	470 x 315 x 250 mm
Weight	4.5 kg
Tabletop Size	500 x 405 mm

Applicable Lighting

Placido's LED lighting	White LED
Fluorescein LED lighting	Blue LED 460 nm
Pupillometry	LED lightingIR LED 875 nm

Notes

Accuracy and repeatability error	Class "A" as per "ISO19980:2005 (E)
Power supply	24V DC external power supply unit
Input power supply unit	90-264 V AC: - 47/63 Hz Max 0.9 A OUTPUT: 24 V DC - 2 A
Power frequency	(50/60Hz) magnetic field IEC 61000-4-8
Power cable	Four-core cable conductors
Computer connection	USB3 Type A cable






Retinal Camera

COBRA PLUS



Retinal Camera COBRA PLUS

1. Head Rest
2. Optical Lens
3. Instrument
4. Fixation Point
5. Chin Rest
6. Chin Rest Hand Held
7. Joystick
8. Guide Guards
9. Phoenix Software (Included)
10. Two Device Table Top (Optional)
11. LIFT 02 Electrical Stand (Optional)
12. Table Up / Down Control
13. Table Storage Drawer
13. Personal Computer (Optional)



Retinal Camera **COBRA PLUS**

The Cobra + is a non-mydiatic fundus camera that comprises all the functions required for a rapid screening of the status of the retina. This ergonomic design provides a clear and detailed image of the ocular fundus with a field of vision of up to 50°.



Advanced Technology

Cobra+ has 9 internal fixation points that allows to capture panoramic images of the peripheral areas. Cobra+ also can record simultaneously colour and infrared images through a CCD high resolution sensor.



High Precision Measurements

Cobra+ uses a manual acquisition and electronically guided joystick to ensure precise focus control and repeatability measurements for multiple fixation points.



Phoenix

(Software Included)

The Cobra Series uses an advance software platform called “Phoenix” which allows the patients data to be saved for future review and analysis.

Multiple wave-length images can be displayed on one screen such as: MGD analysis, choroidal, vascular, nerve fiber, infrared and red-free images.




Data Communication

USB 3.0 connection between the device and the PC enables a fast, easy transfer of the images, and saved in the Phoenix patient management software.




AVR evaluation module (optional)




The AVR tool measures the ratio between the branch arteriolar-venous diameter. A low ratio between the dimension of the vessels, may be predictive of cardiovascular problems in adult patients

Cup to disk measurement




The measurement of the Cup to Disk ratio is easily achieved using the built in measurement tools that are available in the Phoenix software platform for the detection of glaucomatous disease.

Mosaic function



The COBRA+ can capture multiple images up to 50° field of view) which can be combined together in order to create a panoramic image of the peripheral retina.

AVR evaluation module (optional)




Multiple wave-length images can be displayed on one screen: the original image, infrared image red-free image; as well the choroidal, vascular and nerve fiber images.

Integration tool with ERG test



COBRA+ retinal fundus images can be combined with the multi focal ERG test, performed with the RETIMAX device. This new module provides a precise indication of the functionality of every analyzed retinal area (useful for follow-up of Macular Degeneration and degenerative hereditary retinal diseases).

Multiple wave-length images



Multiple wave-length images can be displayed on one screen: the original image, infrared image red-free image; as well the choroidal, vascular and nerve fiber images.

Retinal Camera COBRA HD



Product Specifications

Measurements

Image Resolution	2448 x 2051 (5M Pixel)
Working Distance	20mm
Dimension	420mm (w) x 315mm (l) x 255mm (h)
Shelf Size	380mm (w) x 500mm (l)
Weight	6kg
Base Movement	105mm (w) x 110mm (l) x 30mm (h)
Field of View	50 x 45

Light Source

Auxiliary IR	LED @850nm
White Flash	LED @450-650nm

Notes

Operating Environment	Temperature	-10 °C ~ +35 °C
	Humidity	30% ~ 90% RH
	Atmospheric pressure range	800 hPa ~ 1060 hPa
Storage and Environment Condition	Temperature	-10 °C ~ +55 °C
	Humidity	10% ~ 95%
	Atmospheric pressure range	700 hPa ~ 1060 hPa
Shipping Condition	Temperature	-40 °C ~ +70 °C
	Humidity	10% ~ 95%
	Atmospheric pressure range	500 hPa ~ 1060 hPa
Vibration	10Hz @ 500Hz, 0.5g Shock 30g Duration 6ms,	
	Bump 10g Duration 6ms,	
Power Supply	External power source 24 VCC In: 100-240	
	Vac - 50/60Hz - 0.9-05A - Out: 24Vdc - 40W	
Frequency Range	80MHz - 800MHz	

Additional Accessories

Isolation Transformer	230V / 230V
Motor Driven Telescopic Column	SCHUMO AG, Model TES2 23 / TA0113 X20 400238Z



Retinal Camera

COBRA HD



Features:

- With its ergonomic design Cobra provides a clear and detailed image of the ocular fundus with a field of vision of up to 50 degrees.
- Cobra HD uses the Phoenix software platform allowing patient data to be saved for future review and analysis.
- Cobra HD allows the acquisition of multiple images, to create a panoramic image of the peripheral retinal areas.
- The instrument, with manual acquisition and electronically guided control guarantees high precision and repeatability of the measurements.




Product Introduction

Retinal Camera

COBRA-HD




The **Cobra HD** is a non-mydratic fundus camera that compromises all the functions required for a rapid screening of the status of the retina.



Cobra HD includes a module for the analysis of the Meibomian Glands Using Pheonix soft ware, the glands structure and health can be analysed..




Cobra HD shares the use of the CCD high-resolution sensor (5 megapixel) for the alignment of the patient and the capture of retinal images.



Multiple wave-length images can be displayed on one screen: the choroidal, vascular, nerve fiber, infrared and red-free images,



Cobra HD USB connection between the device and the PC enables a fast and easy transfer of the images and saved in the Phoenix patient management software.




AVR Evaluation Module - The AVR tool measures the relationship between the branch arteriolar-venous diameter.

Retinal Camera

COBRA HD

1. Head Rest
2. Optical Lens
3. Instrument
4. Fixation Point
5. Focus Adjustment Handle
6. Chin Rest
7. Chin Rest Hand Held
8. Joystick with Acquisition Button
9. Guide Guards
10. Personal Computer
11. Electrical Table & Tabletop (Optional)



Retinal Camera COBRA HD



Product Specifications

Measurements

Image Resolution	2448 x 2051 (5M Pixel)
Working Distance	20mm
Dimension	420mm (w) x 315mm (l) x 255mm (h)
Shelf Size	380mm (w) x 500mm (l)
Weight	6kg
Base Movement	105mm (w) x 110mm (l) x 30mm (h)
Field of View	50 x 45

Light Source

Auxiliary IR	LED @850nm
White Flash	LED @450-650nm

Notes

Operating Environment	Temperature	-10 C ~ +35 C
	Humidity	30% ~ 90% RH
	Atmospheric pressure range	800 hPa ~ 1060 hPa
Storage and Environment Condition	Temperature	-10 C ~ +55 C
	Humidity	10% ~ 95%
	Atmospheric pressure range	700 hPa ~ 1060 hPa
Shipping Condition	Temperature	-40 C ~ +70 C
	Humidity	10% ~ 95%
	Atmospheric pressure range	500 hPa ~ 1060 hPa
Vibration	10Hz @ 500Hz, 0.5g Shock 30g Duration 6ms, Bumb 10g Duration 6ms,	
Power Supply	External power source 24 VCC In: 100-240 Vac - 50/60Hz - 0.9-05A - Out: 24Vdc - 40W	
Frequenzy Range	80MHz - 800MHz	

Additional Accessories

Isolation Transformer	230V / 230V
Motor Driven Telescopic Column	SCHUMO AG, Model TES2 23 / TAO113 X20 400238Z;



OSIRIS-T

Aberrometer & Topographer





Features:


- Osiris -T has the ability to measure high order aberrations and topography calculations of the internal component of the wave-front.
- Osiris-T is a total ocular aberrometer, and is indispensable for the correct evaluation of critical patients who have, in addition to traditional low-order defects, even more complex ocular aberrations.
- Osiris-T has a unique design that enables it to measure aberrations with a resolution of 45,000 points (at the maximum pupil diameter), with a wide dynamic.
- Osiris-T is also able to measure the total wave-front in real time with a frame rate of up to 33 images per second: this makes it possible to measure and view changes in power and aberrations while the patient is accomodating.



Product Introduction

Aberrometer + Topographer

OSIRIS-T




Aberrometer & Topographer

Osiris-T has a unique design that enables it to measure aberrations with a resolution of 45,000 points and topographer system based on a 22 ring Placido disk.




Software Phoenix

Osiris-T uses Phoenix software platform allowing patient data to be saved for future review and analysis, shared by all CSO devices.




Dynamic Accomodations

Real-time measurements of the ocular wave-front with customizable exam modes (ramps or square waves) to evaluate patients the ability to focus.



Toric Lens Assistant

Osiris-T aberrometer makes it possible to distinguish whether any astigmatic residue is due to a rotation of the lens or to an incorrect calculation.



Analysis Software for Aberrations

The Phoenix software offers a wide range of analysis options such as refractive error maps and visual stimulations (PSF, MTF and with optotype).



Densitometry


Osiris-T can acquire backlight images without reflections in order for an assessment of cataract and optical media opacity evaluation.

OSIRIS-T

Aberrometer + Topographer

1. Chinrest Module
2. Headrest
3. Capturing Channel
4. Chinrest Support
5. Chinrest Adjuster Knob
6. Handrest
7. Drawer
8. Guide Guard
9. Power In

10. LIFT Electrical Tables
11. Joystick (Capturing Trigger Button)
12. OSIRIS -T instrument
13. Main Supply Cable
14. Switching Adapter (Data Nameplate)



PC: 4 GB RAM - Video Card 1
GB RAM (not shared) resolution
1024 x 768 pixels - USB 3.0 type A
Operating system: Windows XP, Windows
and Windows 10 (32/64 bit).

LIFT 02

Electrical Table (Optional)



Tabletops

Available Colors

OSIRIS-T

Aberrometer & Topographer



Product Specifications

Technical Data

Data Transfer	USB 3.0
Power Supply	External power source 24 VCC In: 100-240Vac 50/60Hz - 0.9-05A Out: 24Vdc - 40W
Power net cable	with plug C14
Dimensions (HxWxD)	515 x 315 x 255mm
Weight	6.9Kg
Chin rest movement	70mm \pm 1mm
Minimum height of the chin cup from the table	24cm
Base Movement (xyz)	105 x 110 x 30mm
Working distance	74mm

Light Sources

Aberrometer	Led @850nm
Fixation	Led @450-650nm
Placido	Led @635nm
Pupillometry and auxiliary	Led @780nm

Aberrometry

Points measured at maximum pupil	45000
Spatial resolution	41 μ m
Pupil size range	2-9mm
Dioptric range	Sph from -25D to +15D; Cyl up to 10D
Repeatability	0.05D on test eyes

Topography

Placido rings	22
Measured points	5632
Topographic covering (in 43D)	10mm
Dioptric measurement range	from 1D to 100D
Measurement accuracy	Class A according to UNI EN ISO 19980-2012
Compatibility with standard	DICOM v3



OSIRIS

MADE IN ITALY
Aberrometer

LIFT 02

Electrical Table (Optional)




Tabletops

Available Colors



Features:

- Osiris has the ability to measure high order aberrations as well as standard refraction has become the new standard of care for your patients.
- Osiris, is a total ocular aberrometer, and is indispensable for the correct evaluation of critical patients who have, in addition to traditional low-order defects, even more complex ocular aberrations.
- Osiris has a unique design that enables it to measure aberrations with a resolution of 45,000 points (at the maximum pupil diameter), with a wide dynamic.
- Osiris is also able to measure the total wave-front in real time with a frame rate of up to 33 images per second: this makes it possible to measure and view changes in power and aberrations while the patient is accomodating.



Aberrometer

OSIRIS

Product Introduction




Ocular Aberrometer

Osiris has a unique design that enables it to measure aberrations with a resolution of 45,000 points (at the maximum pupil diameter), with a wide dynamic.




Software Phoenix

Osiris uses Phoenix software platform allowing patient data to be saved for future review and analysis, shared by all CSO devices.




Dynamic Accomodations

Real-time measurements of the ocular wavefront with customizable exam modes (ramps or square waves) to evaluate patients the ability to focus.



Toric Lens Assistant

Osiris aberrotation makes it possible to distinguish whether any astigmatic residue is due to a rotation of the lens or to an incorrect calculation.



Analysis Software for Aberrations

The Phoenix software offers a wide range of analysis options such as refractive error maps and visual stimulations (PSF, MTF and with optotype).




Densitometry

Osiris can acquire backlight images without refractions in order for an assessment of cataract and optical media opacity evaluation.

OSIRIS *Topographer*

1. Chinrest Module
2. Headrest
3. Shaped Tabletop
4. LED Shooting Lens
5. Main Supply Cable
6. Switching Adapter
7. Data Nameplate
8. Chinrest
9. Handrest
10. Joystick with Capturing Trigger Button
11. Instrument Power Supply Cable
12. Guide guard



OSIRIS

Aberrometer



Product Specifications

Measurements

Operation Distance	79 mm from corneal vertex
Number of Rings	24
Number of Measuring Points	5632
Number of Points Analysed	Over 100000
Diameter of the corneal area covered	0.4 to over 9.6 mm of Diameter
Dioptres Measuring Arc	1 to 100 D
Size (HxWxD) mm	510 x 313 x 280 mm
Weight	6.5 kg

Applicable Lighting

Placido's LED lighting	White LED
Fluorescein LED lighting	Blue LED 460 nm
Pupillometry	LED lighting IR LED 850 nm

Notes

Accuracy & repeatability error	Class "A" as per "ISO19980:2005 (E)
Power supply	24V DC external power supply unit
Input power supply unit	90-264 V AC: - 47/63 Hz Max 0.9 A OUTPUT: 24 V DC - 2 A
Power frequency	(50/60Hz) magnetic field IEC 61000-4-8
Power cable	Four-core cable conductors
Computer connection	USB3 Type A cable



Clinical Slit Lamp
MADE IN ITALY

ZOOM
Magnification

SL-9900





Magnification Options

The SL-9900 Series microscopes has multiple continuous enlargements, giving a clear, brilliant and well contrasted images thanks to the multi-layer anti reflection treatment.



	2X	3X	5X	ZOOM
Binocular Head	Converging Stereoscopic	Galilean magnification changer with converging eyepieces		
Magnification	10x, 16x	10x, 16x, 25x	6x, 10x, 16x, 25x, 40x	6x - 33x Continuous
Light Source	LED	LED	LED	LED
Built-in Yellow Filter	No	Yes	Yes	Yes
Digital Imaging Ready	No	Yes	Yes	Yes
Optional Eyepieces	16x, 26x	No	No	No
Weight (Lbs/Kg)	49 / 22	53 / 24	45 / 20	50 / 23

⊕ ⊖ The SL-9800 & SL9900 slit lamps has a rotation of $\pm 90^\circ$ continuous with a Tabo system and actual visual field of 32-6.2.



Precise Adjusters

Viewlights slit lamps are built with high quality optics and very precise mechanical parts, feasible for all adjustable slit positions designed for comfort and a new experience while capturing and recording the patients eyes.




Maneuver Cross Slide Joystick

The SL-9800 & SL-9900 series are built with a precise, smooth and high quality joystick. It provides lateral, longitudinal and vertical movements during all patients examinations.



Filter Insertion Control + Lever

The optical filters devices controllers can selectively transform light of different wavelengths used in many applications Available Filters: Blue, Red, Gray, Green (Red free).



Lightning Card Red LED Reflex


SL-9800 & SL9900 Series slit lamps are built with a tilting device where the diagonal light sources can be vertically projected up to 20° with 5° intervals. This shows to be very helpful in the horizontal optical observation gonioscopy & posterior eye examination.





Illumination System

SL-9800 & SL-9900 series slit lamp offers an LED illumination lamp designed for accurate and detailed observations. The life of the LED is longer than 50,000 working hours.



High Quality Optics

SL-9800 & SL-9900 Series slit lamps provides high quality optics and excellent image quality for all ophthalmology demands. Our slit lamps features one stereoscopic microscope with 2 magnifications, or one Galilean system with 3, 5 magnifications, or with progressive zoom, complete with screw-out eyepieces.

LIFT 02

Electrical Table (Optional)



Tabletops

Available Colors






Phoenix Software

The SL-9900 Series slit lamps includes a Phoenix software technology which manages to collect data, measurements, captures high quality photos and video for all conducted examinations. Also features an image comparison, drawing tools, user calibration sources, network database, and language converters.



**Viewlight Monitor
Optional**



Model: 100710805


Model: 100710807

LIFT Electrical Tables (LIFT-01 & LIFT-02)

CSO sturdy and secured instrument table features an adjustable height sequence designed for 1-2 optical instrument configurations.



Accessory Drawer



Model: 100710803

Table Tops (Not Included)

All CSO surfaces of these tables are designed to be very rigid with minimum deflection so that the alignment of optical elements remains stable over time



Model: 100710801



HR - Digital Video Camera

The new HR digital video camera has been particularly designed with hardware, firmware and optimized for all ophthalmologic purposes.

Optional - [Accessory](#)




USB CONNECTION

HR camera is based on a high performances CCD sensor, an excellent color rendering and a integrated Phoenix software.



Canon Cameras

All our hightech slit lamps offers an optional Canon EOS Rebel Series cameras deisgned to capture photos and videos with its high resolution megapixel sensors and precise focus.



Optional Adaptors


Viewlight USA slit lamp are designed to incorporate optional accessories such as the beam splitter and camera adaptor so that it can be attached with the stylish Canon professional camera.

Digital Adaptors [VL-Z8](#)



Tonomter Z-800

The Z-800 tonometer provides extreme measurement accuracy. All values are directly readable on the instrument and without any standardization or calibration difficulties.




Tonomter A-900

The A-900 is an accessory to any slit lamp model, used for measuring ocular pressure, amount of force needed to temporarily flatten part of your cornea and to prevent risk of developing glaucoma.



A-900

Tonometer




Product Specifications

Measurements

Pressure Surface	7mm
Measurement Drum Surface	3.06mm (7,354mm area)
Conversations Pressure	mmHG to Kpa
Angle Light Source	60°
Tonomoter Scale	6 Calibrations
Measurement Force	Generated by the Spring
Measurement Range	0 / 80 mmHG (0/10,64kPA)
Weight	0.48 kg (without accessories)

Z-800

Tonometer



Product Specifications

Measurements

Pressure Surface	7mm
Measurement Drum Surface	3.06mm (7,354mm area)
Conversations Pressure	mmHG to Kpa
Angle Light Source	60°
Tonomoter Scale	6 Calibrations
Measurement Force	Generated by the Spring
Measurement Range	0 / 80 mmHG (0/10,64kPA)
Weight	0.48 kg (without accessories)

Slit Lamp SL-9900 Series

SL-9900 Series are designed slit lamps for ophthalmologists and optometrists for specific diagnostic procedures such as the biomicroscopic examination of the eye, ocular fundus, and posterior vitreous body.

1. Light Bulb / LED Cover
2. Slit Height Value Index
3. Extractable Eyepieces
4. Filter Selection Knob
5. Slit Height Adjuster / Tuner
6. Microscope Splitter Knob
7. Microscope Locking Knob
8. Fluorescein Filters Insertion Rod

9. Projector Positioning Scale
10. Slit Width Adjustment Knob
11. Brightness Control Knob

12. Joystick (lateral, longitudinal & vertical)

13. Shaped Table Top (Optional)

14. Fixation Point Power Supply

15. Chinrest Module

16. Raduated Scale

17. Head Rest

18. Magnification Tuner

19. Slit Lamp Mirror

20. Fixation Point

21. Microscope Arm Fixing Knob


22. Chin Rest

23. Chinrest Height Adjusting Ring Nut

24. Slit Width Adjustment Knob

25. Patient's Handle

26. Wheel Shields




SL-9800 & SL9900

Slit Lamp Series

Specifications

Stereoscopic Microscope

Microscope	2x	3x	5x	Zoom
Type	Convergent	Galleian convergent with magnification change syste		Galleian convergent w/ variable magnification
Eyepiece Convergence Angle	13°		6°	
Eyepieces	10x	12,5x	12,5x	12,5x
Eyepiece adjustment	±8 D.			
Magnifications	10x 16x	10x 16x 25x	6x 10x 16x 25x 40x	7X 30X
Field of View		8,5x 14,8x 25,6x	5,6x 8,5x 14,8x 25,6x 39,3x	
Interpupillary distance	18,5mm 12mm	26mm to 8,5mm	41mm to 5,7mm	30mm to 7,4mm
Barrier filter	51,5mm to 87mm	50mm to 80mm		
Barrier filter	Yellow			



Technical Data

Slit Projection	1X
Slit Width (setting continuous)	0 – 12 mm
Slit length (setting continuous)	1– 12 mm
Slit Length (max)	12 mm
Apertures	12, 9, 5, 1, 0.2 mm
Filters	Blue, Red, Gray, Green (Red free)
Light Diffuser	Light Diffuser
Background Light	Only with Digital Vision HR
Slit Rotation	± 90° continuous (Tabo system)
Slit Angle	Variable 0° 5° 10° 15° 20°
Rotation - Slit projection	±90°, angular scale, ref. on 0° and ±10°
Working distance	80 mm
Joystick Push Button	Only with Digital Vision HR
Left / Right Detection	Only with Digital Vision HR
Voltage	15V DC 1A
Light Source	White LED
Brightness	248000 LUX conti nuous adjustment
Dimensions (HxWxD)	675 x 313 x 335mm
Weight	7 kg

System Requirements (Version Digital HR)

PC:	4 GB RAM
Video Card:	1 GB RAM
Resolution:	1024 x 768 pixels
USB	3.0 Type A
Operati ng System	Windows XP, Windows 7 and Windows 10 (32/64 bit).



SLIT LAMP SL-9900



Elite



LIFT 02

Electrical Table (Optional)



Tabletops

Available Colors





SLIT LAMPS

The new SL-9900 ELITE slit lamps, works with an innovative LED illumination system. Offers a new beam splitter equipped with the new high resolution digital camera, and it has been designed for whom wants high performances and smart design.



5^x

Magnification Tuners

The SL-9900 Elite slit lamps have multiple continuous enlargements, giving a clear, brilliant and a well contrasted images thanks to the multi-layer antireflection treatment.



Galilean magnification changer with converging eyepieces

Magnification	6x, 10x, 16x, 25x, 40x
Light Source	LED
Built-in Yellow Filter	Yes
Digital Imaging Ready	Yes
Optional Eyepieces	No
Weight (Lbs/Kg)	45 / 20



Maneuver Joystick

The SL-9900 Elite are built with a precise, smooth and high quality joystick. It provides lateral, longitudinal and vertical movements during all patients examinations. Enables you to optically change the illumination field up by a single click actions; ideal for contact lens fitting observation.



Precise Slit Lamp Adjusters

The L-9900 Elite slit lamps are built with high quality optics and precision mechanical parts feasible for all adjustable slit positions designed for comfort and a new experience while capturing & recording all examinations.





Illumination System

The advanced SL-9900 Elite offers a new LED illumination source therefore allowing the examination with very comfortable and bright experience while taking clear pictures with each magnification device.

1. LED Power Illuminator
2. LED Background Illumination (Built-in)
3. Lighting Card Red LED Reflex

Lightning Card Red LED Reflex

SL-9900 Elite slit lamps are built with a tilting device where the diagonal light sources can be vertically projected up to 20° with 5° intervals. This shows to be very a helpful horizontal optical observation gonioscopy & posterior eye examination.



Filter Insertion Control + Lever

The slit lamp SL-9900 Elite offers an optical filters devices controllers can selectively transform light of different wavelengths used in many applications available filters: Blue, Red, Gray, Green (Red free).




Illumination System

The SL-9900 Elite offers an LED illumination lamp designed for accurate and detailed observations. The life of the LED is longer than 50,000 working hours and its color temperature is constant in each situation.



High Quality Optics

L-0990 Elite slit lamps provides excellent high optics and excellent image quality for all ophthalmology demands. These converging microscope optics provides comfort to the user when using the slit lamp. All the microscopes have a yellow filter to improve the image quality during fluorescein tests.




Specifications:



Slit Lamp SL-9900 ELITE

1. Light Bulb / LED Cover
2. Slit Height Value Index
3. Filter Selection Knob
4. Slit Height Adjuster / Tuner
5. Microscope Splitter Knob
6. Extractable Eyepieces
7. Digital Video Camera
8. Joystick (lateral, longitudinal & vertical)
9. Brightness Control Knob
10. Orthogonally Moving Base
11. Raduated Scale
12. Magnification Tuner
13. LED Background Illumination
14. Magnification Tuner
15. Lighting Card Red LED Reflex

16. Microscope Locking Knob
17. Projector Positioning Scale
18. Eye Positioning Reference Index
19. Head Rest
20. Fixation Point LED Illuminator
21. Chin-Rest
22. Tonometer Plate
23. Chin-rest Height Adjuster
24. Microscope Arm Fixing Knob
25. Projector Arm Fixing Knob
26. Patient's Handle
27. Wheel Shields
28. LED Light Cable
29. Power Supply Cable
30. Device Base Lock





HR Digital Video Camera

The new HR digital video camera has been particularly designed with hardware, firmware and optimized for all ophthalmologic purposes. HR camera is based on a high performances CCD sensor, an excellent color rendering and a integrated **Phoenix Software**.




Canon Cameras

SL-9900 Series slit lamps are designed to incorporate optional accessories such as the beam splitter and camera adaptor so that it can be attached with the stylish Canon professional camera.

1. Digital Camera Beam Splitter (Optional)



Phoenix
(Software Included)



HR-Elite Digital Camera



Product Specifications

Image Sensor Type 1/1.8" progressive scan color CCD

Picture Size (LIVE) Up to 1624 (h) x 1232 (v) @8fps

Cell Size 4.40 μm x 4.40 μm

Resolution depth 14 bit

Digital Interface IEEE 1394a; DCAM V1.31

Transfer Rate 400 Mb/s

Frame Rates 15 fps

Video Modes 1600 x 1200, 800 x 600, ROI

Power Consumption Less than 3 W

Minimum PC Requirements

Processor i5 or newer

RAM 8GB

Hard Drive 1 Tera Byte

Connection USB 3.0

Resolution 1280x1024 pixels

Operating System Microsoft Windows XP , Microsoft Windows 7,
Microsoft Windows 10

Networking TCP/IP protocol for network



MADE IN ITALY

Slit Lamp SL-9900

Elite

Stereoscopic Microscope

Microscope 5x	
Type	Galilean convergent with magnification change system
Eyeiece Convergence Angle	6°
Eyeieces	12,5x
Eyeiece adjustment	±8 D.
Magnifications	6x 10x 16x 25x 40x
Field of View	8,5x 14,8x 25,6x
Interpupillary distance	41mm to 5,7mm
Barrier filter	50mm to 80mm
Barrier filter	Yellow



Technical Data

Slit Projection	1X
Slit Width (setting continuous)	0 – 12 mm
Slit length (setting continuous)	1– 12 mm
Slit Length (max)	12 mm
Apertures	12, 9, 5, 1, 0.2 mm
Filters	Blue, Red, Gray, Green (Red free)
Light Diffuser	Light Diffuser
Background Light	Only with Digital Vision HR
Slit Rotation	± 90° continuous (Tabo system)
Slit Angle	Variable 0° 5° 10° 15° 20°
Rotation - Slit projection	±90°, angular scale, ref. on 0° and ±10°
Working distance	80 mm
Joystick Push Button	Only with Digital Vision HR
Left / Right Detection	Only with Digital Vision HR
Voltage	15V DC 1A
Light Source	White LED
Brightness	248000 LUX continuous adjustment
Dimensions (HxWxD)	675 x 313 x 335mm
Weight	7 kg

System Requirements (Version Digital HR)

PC:	4 GB RAM
Video Card:	1 GB RAM
Resolution:	1024 x 768 pixels
USB	3.0 Type A
Operating System	Windows XP, Windows 7 and Windows 10 (32/64 bit).




Phoenix Software

The SL-9900 - Elite slit lamps includes a Phoenix software technology which manages to collect data, measurements, captures high quality photos and video for all conducted examinations. Also features an image comparison, drawing tools, user calibration sources, network database, and language converters.



**Viewlight Monitor
Optional**



Model: 100710807

Model: 100710805

LIFT Electrical Tables (LIFT-01 & LIFT-02)

CSO sturdy and secured instrument table features an adjustable height sequence designed for 1-2 optical instrument configurations.



Accessory Drawer



Model: 100710803

Table Tops (Not Included)

All CSO surfaces of these tables are designed to be very rigid with minimum deflection so that the alignment of optical elements remains stable over time



Model: 100710801

NEW LIFT

01 & 02

Electrical Tables

.....
Table Tops - Not Include



LIFT *ELECTRICAL TABLES*

Overall Features

The new ophthalmic electrical table **LIFT 01 & 02** allows to place one or more devices on its table top, whose elevation can be adjusted. The electric table is composed of an optional table top where, cogged guides are already installed for the device housing.

The “strong” design, obtained by the increasing of weight and dimension, the use of new material (no plastic), gives a modern and essential appearance.

- Power cable connection
- Cogged Wheels
- Scrolling plate
- Keyboard Connection
- Speed: 10mm/s
- Voltage: 100-240V



Electrical Table Elevation




Table top elevation can be adjusted by using the keypad which operates on the elevation telescopic column.

Smooth Transitions


Multivoltage elevation column from 100Vac to 240Vac (50Hz o 60Hz) and offers a power supply cable (German and USA plug).

NEW

LIFT - 01 Electrical Table

Standard base is made by die-cast aluminum. Its weight is 6,5kg but can be increased using ballasts to fix under the base (+7kg). It is painted by anti-scratch coating. It can support feet or wheels.

- Power Cable Connection
- Cogged Wheels
- Max Thrust Load: 1400N
- Speed: 10mm/s
- Voltage: 100-240V
- Absorption Max: 1,5A



Multi Voltage - Elevation

Multivoltage elevation column from 100Vac to 240Vac (50Hz o 60Hz) with wheels and feet power supply cable (German and USA plug) Keypad for up/down control with cable.



LIFT - 01


Table Tops

STANDARD DESIGN - Optional

Standard base is made by die-cast aluminum. Its weight is 6,5kg but can be increased using ballasts to fix under the base (+7kg). Painted by anti-scratch coating. It can support feet or wheels.

- Made out of Wood
- Tools Compartment
- Cogged Wheels
- Max thrust load: 1400N



Model: 100710801



40" x 50"



Model: 100710803



Recommended Devices:

1. EndoThelial Microscope - Perseus
2. Topographer & Tomographer - SIRIUS
3. Aberrometer & Topographer - OSIRIS-T
4. Digital Slit Lamp - ELITE





NEW

LIFT - 02 Electrical Table

Lift 02 - is features a adjustable configuration, designed with metal sheets to improve a stable support for 1 or 2 instruments display. Its rugged construction and simplicity allows you peace of mind in knowing your instruments are secured on a stable platform.

1. Stability
2. Hidden cables
3. Easy to assembly
4. Anti-Scratch Coating
5. Heavy (23Kg)




LIFT - 02 Table Tops

MODEL: **100710805**

HANDICAP DESIGN - Optional

The LIFT - 02 table tops are all composed with two different side: Top Side is assembled on the bottom side. It is painted by anti-scratch coating, and generally in white color.

Bottom Side is assembled on the column. It can be paint in different colors. It is wider and stronger respect the other.



MODEL: **100710807**



Recommended - Devices

1. Slit Lamps SL-9900 Series
2. Topographer - Anatares
3. MS-39 Anterior Segment OCT
4. Retinal Camera - COBRA




POLARIS

Tear Film Analysis




POLARIS


Tear Film Analysis



Tear Film Examination



Pupillography



Features:

- Polaris allows you to assess stability and regularity of the tear film, using non-invasive break up time measurement (NIBUT)
- Display and assess tear film lipid patterns
- Analyses the interference between the white light and fringes
- The 20 tilted beam splitter is an available accessory, suitable for all slit lamps models
- Evaluates tear meniscus

Descriptions



POLARIS

Tear Film Analysis



Product Specifications

Measurements

Dimensions	22mm (w) x 217mm (l) x 74.5mm (h)
Weight	300 g
Manual	Alignment
LED Lighting Source	White LED
Computer Connection	USB Cable


Notes

Power Supply Voltage	5V - 100mA
Operating Environment	Temperature +10 C ~ +35 C Humidity 30% ~ 90% RH Atmospheric pressure range 800 hPa ~ 1060 hPa
Storage and Environment Condition	Temperature -10 C ~ +55 C Humidity 10% ~ 95% Atmospheric pressure range 700 hPa ~ 1060 hPa
Shipping Condition	Temperature -40 C ~ +70 C Humidity 10% ~ 95% Atmospheric pressure range 500 hPa ~ 1060 hPa
Vibration	10Hz @ 500Hz, 0.5g Shock 30g Duration 6ms, Bumb 10g Duration 6ms,
Power Supply	External power source 24 VCC In: 100-240 Vac - 50/60Hz - 0.9-05A - Out: 24Vdc - 40W
Frequenzy Range	80MHz - 800MHz



PERSEUS

Endothelial Microscope



NOT SOLD IN
US



Endothelial Microscope PERSEUS

What does PERSEUS do?

PERSEUS is a fully automated non-contacted instrument providing analysis of the corneal endothelial while producing high definition quality and well contrasted images.

What type of diagnoses PERSEUS can detect?

The endothelial microscopy is essential in the diagnosis of many corneal dystrophy, degenerative diseases of pre and post operative assessment of cataract surgery and corneal transplants.

What other medical features it displays?

PERSEUS enables mapping and measurement of the endothelial cells and the acquisition of a series of parameters to evaluate the cornea's health status. See below for all relevant data:


- Cell density and average area
- Automatic focus of the endothelial layer
- Automatic calculation of cells centres and extensive statistical analysis based on the collected data
- Variance coefficient
- Average median error
- Cell size occurrence histogram



Endothelial Microscope PERSEUS

Mosaic Functions

PERSEUS capture multiple images while viewing on one screen 1 central and 6 peripheral



Quick Acquisition of Images

The digital CCD camera enables PERSEUS to automatically track and focus on the patients eye, taking high quality images. Simply touch the LCD screen and PERSEUS will simply do the rest.

- Photographic method : Contactless
- Photographic field: 0.54 mm x 0.27 mm
- Measurement accuracy: $\pm 10 \mu\text{m}$

Comparison Function

PERSEUS application software allows the comparison of all patients so that disease progression can be properly monitored.

Analysis and Detection

- Allows for non-invasive non-contact examination.
- No risk of transmission of infectious diseases
- Pain-free examination not requiring local anaesthesia.
- Non-invasive exam of the endothelial tissue

LIFT 02

Electrical Table - Not Included



Optional Tabletops

Available Colors



PHOENIX - Patient Management Software

PERSEUS application software can be connected to the Phoenix patient management platform, allowing patient data to be saved for future review and analysis.






Non-invasive examination of the Endothelial Layer

PHOENIX software measures and displays a color coded map of the endothelial layer along with a set of indices based on the shape and size of cells.

Evaluation and Manual Editing

PERSEUS enables single acquisition to automatically measure and counts up to 400 cells. It can also calculate the cell density, Pleomorphism, Poimegatism and Pachymetic data values.

- Hexagonal deviation (percentage of hexagonal cells)
- Cell shape occurrence histogram.
- Shape factor



LIFT 01

Electrical Table - Not Included



Optional Tabletops

Available Colors



Endothelial Microscope

PERSEUS

1. Headrest
2. Optical Unit
3. Chinrest
4. 2 USB Ports
5. Ethernet Port
6. Phoenix Software (Included)
7. 10.4" LED Touch Screen
8. Table Top (Optional)
9. Start Button





Measurements

Picture Size	0.54mm x 0.27mm
Measurement Accuracy	±10 µm
LCD screen	10.4"
Dimension	250mm (w) x 448mm (l) x 437mm (h)
Weight	15 kg

Technical Data

Acquisition Mode	Non-contact, Automatic
Camera	CCD Camera
Focusing Light	LED
Magnification	180X
Pachymetry Measurement	from 400 µm - 700 µm step 10 µm increments
Fixation Target	Internal LED

Notes

Power Supply Voltage	100V-120V ac ±10% 230V-240V ac ±10% 50 Hz, 60Hz
Fuses	Mains socket unit: 5x20 mm 2x 1.25AT
Power Absorbed	200 VA
Power Cable	Three-core cable (with protective earth), conductors minimum cross-section 1 mm ²
External Plugs:	Ethernet, 2 USB
Max Power Absorption:	100 VA