






Mizar















While the biomicroscopic examination, the **Phoenix** application software, allows the display of the smallest details and the acquisition of high-resolution images, videos, and plenty of storage data management for further review.



✓ Image and Video Processing
✓ Patient Management Center
✓ User Friendlt Interface
✓ Export and Print Functions
✓ Image and Drawing Tools



MIZAR DIGITAL CAMERA

Resolution	5Мрх
Pixel resolution	2448x2048
Camera interface	USB 3.0
Frame Rate	35 fps
Image Sensor Rate (Size)	2/3"

CPU: I3 or higher (suggested I5), RAM: 4 GB or higher (suggested 8 GB), OS: Win10x64, Win8x64, Win7x32, Win7x64, GPU: 1 Gbyte not shared, RESOLUTION: 1280x960 or higher

