

Opera

OPHTHALMIC UNIT



FEATURES

- Unique styling and functionality
- Manual pneumatic assisted tilt recline examination chair
- Release bar located on chair back with chair base controlled by footswitch.
- Stationary pillow headrest standard.
- CSA, C/US Certificate of Compliance

About

Luxurious Oversized Headrest

Large, padded headrest is fixed in place to keep the patient comfortable in the reclined position.

Comfortable and convenient armrests

Convenient ingress and egress from either side of the chair as upholstered arms rotate upward. The armrests are upholstered and padded for patient comfort.

Maximum Recline Control to 40 degrees

A recline bar, which is easily grasped from both sides of the chair, controls the assisted recline of the chair. The chair reclines to the maximum angle of 40 degrees.

